medivon.

Gun massager

Professional percussion therapy

PL | EN | DE | FR | ES | PT | IT | CZ | SK
SL | EE | LT | LV | HU | RO | RU | NB | FI
EL | HR | SV | DA | NL | TR | BG | AR

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1. Device operation

The device is very easy to use - just one press of a button is enough to start a relaxing full body massage that will allow you to relax after a hard day at work or effectively massage tense muscles after training.

After starting with the power button on the back of the device, you can change the massage modes cyclically, gradually increasing the frequency and power of the vibrations using the + and - buttons

To turn off the device, slide the power button back to the OFF position

Before use, it is recommended to fully charge the device with the charger. Only the charger included in the set should be used for this purpose. A fully charged battery will be indicated by the backlight. It is recommended not to fully discharge the device to extend battery life. It is forbidden to use the massager while charging the battery. It may cause damage or fire.

Note: we replace the massage heads only when the device is turned off. It is enough to remove them from the device by pulling (they do not have a screw thread).

2. Massage and settings

Start by choosing a massage head that is suitable for you and your needs. Make sure it is firmly attached and not damaged. Then start the device by pressing the POWER button.

The strength of the massage is regulated by the user. Depending on the force with which we press the working massager to the body, we will get different final effects. For the first massage, it is recommended to use the massager gently.

Using the M button, we start the heating mode that is adapted to the heating head (the head must be pressed to the very end of the mounting bracket). In this mode, use the + and - buttons to select the temperature of the head in the range of 35 to 60 degrees. The cooling head must be pre-cooled for about an hour in the freezer.

3. Dangers

1. Do not use spare parts that are not approved by the manufacturer.
2. Unprofessional service should not disassemble and replace parts of the leg massager without permission.
3. If the power cord is damaged, it must be replaced by appropriate professional personnel to avoid a hazard.
4. Do not use or store the product in a place with high humidity or high temperature.
5. After completing the massage, disconnect the power cord.
6. Do not use the product for more than 15 minutes at a time.
7. If any side effects occur, please stop using the massager immediately.
8. If the temperature is too high or the product starts to smell bad or smells burnt, stop using it immediately and disconnect the power.
9. Do not carry the product by pulling on the power cord to prevent the power cord from loosening.
10. This product is not intended for people who have no physical, sensory or mental abilities or have no experience and knowledge in this field (including children).
11. People insensitive to heat should use the heating function carefully - preferably under the supervision of a caregiver.

4. Warnings

If you have any of the following conditions, consult your healthcare practitioner before use:

- Pregnant women, children or the elderly
- People who suffer from skin infections or have reddened and burned skin
- who use a pacemaker or similar implanted medical device;
- Patients with heart disease or malignant tumors, etc.
- People who are treated or who experience discomfort
- People suffering from osteoporosis
- People with high blood pressure or abnormal blood pressure

5. Cleaning and maintenance

Turn off the product and unplug the power cord before cleaning. Wipe the surface of the product with a dry, soft cloth. Do not clean the product and accessories with harsh detergents such as benzene, thinners, and store in a cool, dry place when not in use. Also avoid exposure to direct sunlight or humid environment, otherwise it may cause color fading or product damage.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after August 13, 2005.

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The manufacturer hereby declares that the Gun Bandit device complies with Directive 2014/35 / EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, Directive 2014/30 / EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65 / EU of the European Parliament and of the Council of 8 June 2011 on restriction of the use of certain hazardous substances in electrical and electronic equipment. The full text of the declaration of conformity is available at the internet address: www.medivon.pl/deklaracje

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska
Model: Gun Bandit