medivon.

User manual

Massage collar

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1. Device operation

1. Open the valve cap. Press the air filling valve continuously until the air pressure meets the user's personal demand, then close the valve cap.
2. Place the device around your neck and stabilize it properly.
3. Press and hold the power button for 2 seconds, the product will turn on. Press and hold for 2 seconds and the power will turn off automatically.
4. To change the intensity level, briefly press the power button. Low intensity 45°C - orange light, high intensity 50°C - red light.
5. When the product is low battery, the indicator light will flash blue for 5 seconds and then turn off automatically.
6. Blue light flashes while charging. Once fully charged, the blue light remains on continuously.

Before using the device for the first time, fully charge the battery. The first charge may take longer than usual. Using the massager while charging is prohibited and may damage the device or battery. It is recommended not to fully discharge the massager to extend the battery life.

2. Warnings

If you have any of the following conditions, consult your doctor before use:

- pregnant women, children or elderly people,
- people suffering from skin infections or skin redness and burns,
- people who use a pacemaker or similar implanted medical device,
- patients with heart disease or malignant tumors, etc.
- people who feel discomfort,
- people suffering from osteoporosis,
- people with high or abnormal blood pressure.

3. Dangers

1. Do not use replacement parts that are not approved by the manufacturer.
2. Non-professional service should not disassemble and replace parts of the device without permission.
3. If the power cord is damaged, it must be replaced by appropriate professional personnel to avoid danger.
4. Do not use or store the device in a place with high humidity.
5. Make sure that the massager is not exposed to direct sunlight or placed near radiators or other heat sources.
6. Do not use the product on a wooden surface, the zipper may scratch the surface.
7. If any side effects occur, stop using the device immediately.
8. If the temperature is too high or the product begins to smell bad or emit a burning smell, stop using it immediately and unplug the power.
9. This appliance is not intended for use by persons who lack physical, sensory or mental capabilities or lack of experience and knowledge (including children).
10. People who are not sensitive to heat should use the heating function carefully - preferably under the supervision of a caregiver.

4. Cleaning and care

Turn off the product before cleaning. Wipe the surface of the product with a dry, soft cloth. Do not clean the product and accessories with corrosive detergents such as benzene or thinners. Store in a cool, dry place when not in use. Also avoid exposure to direct sunlight or humid environments, otherwise it may cause colors to fade or damage the product.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after 13.08.2005.

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The manufacturer hereby declares that the Medivon CosyGo device complies with Directive 2014/35/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65/EU of the European Parliament and of the Council of June 8, 2011 on restrictions on the use of certain hazardous substances in electrical and electronic equipment. The full content of the declaration of conformity is available at the following internet address: www.medivon.pl/deklaracje

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Importer: Armare S.A. | ul. Grzybowska 87, Warszawa, 00-844, Polska