medivon®

User manual

EYE MASSAGER

PL | EN | DE | FR | ES | PT | IT | CZ | SK
SL | EE | LT | LV | HU | RO | RU | NB | FI
EL | HR | SV | DA | NL | TR | BG | AR

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Weight: around 750 g

Normal working conditions:

Operating environmental temperature: +5°C~+40°C

Relative humidity: ≤80%

Atmospheric pressure: 86kPa~106kPa

Transportation and storage conditions:

Storage environmental temperature: -20°C~+55°C

Relative humidity: ≤93%

Atmospheric pressure: 50kPa~106kPa

2. Device operation

A short press of the POWER(3) button turns on the device. If the battery LED flashes after startup, please charge the device using the included USB cable. It is recommended to fully charge the battery before using the device for the first time. After putting the massager on your head, make sure it is positioned comfortably, you can adjust the settings using the elastic band (1)(2).

Device modes are controlled by repeatedly pressing the power button (3). Subsequent functionalities that are currently active will be highlighted on the screen using LEDs (6). The device has 5 individually prepared cyclical operating modes.

The eye massager has a Bluetooth function that is enabled by default when the device is turned on. Just connect with your mobile device and search for "Horizon Pro" from the list.

A short press of the music button (4) starts the next song from the list, and a longer press (4) turns off the music mode.

To turn off the massager after relaxation, hold down the POWER(3) button again.

3. Warnings and contraindications

Before using the product, please read and understand the operating instructions carefully.

Please clean the product's eye shield before use.

Children and people unable to express their needs should not use the product.

People with sensitive skin or insensitive to heat should not use the product.

Pregnant women, people suffering from heart disease, circulatory disorders, malignant tumors, cerebrovascular diseases and other acute diseases, as well as people under the care of a doctor, should not use the product.

If you experience any discomfort or abnormal skin reaction while using the product, stop using it immediately and consult a doctor.

People with eye problems should avoid using the product or use it under medical supervision. People who have undergone eye surgery and patients suffering from glaucoma, cataracts, retinal detachment and other eye diseases should not use the product.

People wearing contact lenses should remove them before using the product. Do not use the product in places with high temperature, flammable, humid or places with electromagnetic radiation.

Do not disassemble the product without permission; Please make sure the product is free from impact and water to avoid failure.

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The manufacturer declares under his sole responsibility that the Horizon Pro device complies with Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, as well as Directive 2011 of the European Parliament and of the Council /65/EU of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to the Directive of the European Parliament and 2011/65/EU as regards the list of restricted substances. The full text of the declaration of conformity is available at the internet address: www.medivon.pl/deklaracje

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after August 13, 2005.

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Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska