

SURESIGN COVID-19 Rapid Test

- COVID-19 Antigen Rapid Test (Oral Fluid)
- Charting the course of COVID-19
- Why quick, easy saliva-based antigen tests will play an important role.
- Accuracy: 98.9%
- Detects virus as low as 100 pg/ml

What are the advantages of using SURESIGN Antigen Rapid Test (oral fluid)?

Saliva enables non-invasive frequent testing and reduces the need for trained healthcare professionals during collection. SURESIGN COVID-19 Antigen Rapid Test has numerous advantages: the saliva specimen is easily self-collected by the patient with almost no discomfort (allows for effective testing of children, elderly and fragile patients), it does not require specialized health care personnel for its management, and it reduces the risks of exposure for the operator.



Article number: 1875-00550

EAN-number: 5060100937687

Cases: 185x130x80mm – 453g, 20 tests

Cardboard: 550x395x460mm – 14.5kg, 32 cases x 20 = 640 tests

Cardboard EAN-number: 5060100937670

Transport packaging: Pallet 1200x1000x1530mm - 232.6kg, 15 Cardboard = 9600 tests

Country of origin: Northern Ireland

